H-spitalSafety SOLUTIONS

UVC deactivates
E. coli, SARS,
MRSA & COVID-19

Transform Your Pass Thru Window



From This







- Custom Built
- Green Technology
- No Harsh Chemicals
- Safe UVC Disinfection
- Installation Available
- 1-Minute Load
- Reproducible Results
- Reliable Medical Usage since 2009

Why further disinfect hand-washed items?

Hand-washed items are still dirty.

A study published in a peer-reviewed journal,

AAMI/Horizons, "The average cfu's of Bioburden was
reduced from 36.6 cfu's to .3 cfu's after 1 minute of UVC."



Contact Us:



hospitalsafetysolutions.com



844-889-6757



info@hospitalsafetysolutions.com



US Patents# 9,295,741 B2 8,143,596 B2

HespitalSafety SOLUTIONS



How do YOU disinfect "high touch" items?

UVC DISINFECTS in just 60 seconds - 99.999% effective

- Cell Phones
- Tablets/PDA's
- Laptops
- Badges
- Masks

- Eyewear
- Glasses
- Loops
- B.P. Cuffs
- Stethoscopes



At 254 nm, UVC light breaks the molecular bonds inside microbial DNA, thereby inactivating them.



Standard Countertop Unit

- Easily fits three open 15" laptops
- Removable shelves for larger items
- Height: 36", Width: 23", Depth: 29"
- Weight 125 lb.

Contact Us:



hospitalsafetysolutions.com



844-889-6757





US Patents# 9,295,741 B2 8,143,596 B2